



ECL Detection Kit (Cat. No.: S1GNM03j20001)

Description:

ECL Detection Kit is a horseradish peroxidase (HRP) substrate and it is widely utilized for the identification of specific proteins. As shown in Fig. 1, protein samples are initially separated through SDS-PAGE and then transferred onto either nitrocellulose or polyvinylidene fluoride (PVDF) membranes. Following a blocking step, the membrane undergoes probing with a primary antibody. After a thorough washing procedure, the membrane is incubated with a secondary antibody conjugated to horseradish peroxidase (HRP) enzyme. After a final washing step, the membrane is exposed to the ECL Detection Kit to generate a detectable signal.




Fig 1. Principles of ECL Western blotting Protein Detection

Kit Contents:

Contents	S1GNM03j20001 (50 mL x 2)
ECL Detection Kit, Reagent A	50 ml x 1 bottle
ECL Detection Kit, Reagent B	50 ml x 1 bottle

**Storage:**

This product can be shipped at room temperature. For long-term storage, the ECL Detection Kit is stable for 6 months at 4°C when protected from light.

Protocol:**1. Blocking:**

- After protein transfer, remove the membrane from the transfer device.
- Incubate the membrane with blocking buffer to block nonspecific binding sites for 1 hour at room temperature with shaking or overnight at 4°C without shaking.

2. Primary Antibody Incubation:

- Remove blocking solution and add the diluted primary antibody solution.
- Incubate for 1 hour at room temperature with shaking.

3. Washing:

- Wash the membrane three times with TBST wash buffer (TBS with 0.05-0.1% Tween 20), each for five minutes, to remove unbound antibody.

4. Secondary Antibody Incubation:

- Incubate the membrane with diluted secondary antibody solution (HRP conjugated) for 1 hour at room temperature with shaking.

5. Washing:

- Wash the membrane three times with TBST, each for five minutes, to remove unbound antibody. Additional washes can help reduce background.

6. Substrate Preparation:

- Mix equal volumes of Reagent A (black plastic bottle) and Reagent B (white plastic bottle) to prepare the substrate working solution.
- Prepare only enough working solution to cover the membrane (e.g., 6-7ml per 10 X 5 cm membrane).

Note: Use the prepared substrate working solution immediately after mixing. It remains stable for up to 1 hour at room temperature.

7. Substrate Incubation:

- Incubate the membrane for 1–5 minutes at room temperature without shaking. Ensure the working solution covers the membrane.

**8. Membrane Handling:**

- Remove the membrane from the solution, then use an absorbent towel to blot excess liquid, and place it in a plastic sheet protector.

9. Exposure:

- In a dark room with a safe light, position the covered membrane in a film cassette with the protein side facing up.
- Place the film on top of the membrane and expose it for 1 minute.
- Exposure time can be adjusted for optimal results. Note that light emission is most intense immediately after substrate incubation, decreasing significantly within a few minutes.

◎ Revision History ◎

Description	Version	Date
Initial Release	S1GNM03j20001_Protocol_V1	Sep 2023