



First Strand cDNA Synthesis Kit II

(Cat. No.: S2GNM03j30002)

Description:

First Strand cDNA Synthesis Kit II contains recombinant MMLV reverse transcriptase with improved thermostability and reduced RNase H activity. It is an easy-to-use kit for reliable cDNA synthesis, which is able to synthesize the first strand cDNA at 37~50°C.

This kit includes RNase Inhibitor, which effectively inhibits RNase A, RNase B, and RNase C activities, ensuring the integrity of RNA during the reverse transcription process. The product comes with an optimized RT Buffer and an Oligo (dT)/Random Primer Mix, enabling the highly efficient synthesis of short-chain cDNA that is suitable for real-time PCR applications.

Kit Contents:

Contents	S2GNM03j30002 (100 rxns)
RTase/RI Enzyme Mix	100 µl x 1
5X RT Buffer (DTT/dNTPs)	500 µl x 1
Oligo (dT)/Random Primer Mix	100 µl x 1
DEPC-Treated H ₂ O	1 ml x 2

Storage:

The kit is stable for 24 months at -20°C.



Protocol:

1. After thawing, mix and briefly centrifuge the components of the kit, keep it on ice.
2. Add the following reagents into a PCR tube and keep it on ice.

Table 1. Reaction Setup - Denature (Mixture A)	
Components	Volume
RNA template	X μ l (1ng~2 μ g)
Oligo (dT)/Random Primer Mix	1 μ l
DEPC-Treated H ₂ O	to a final volume of 10 μ l

Mix well and heated at 70°C for 5 minutes in advance, then incubated on ice bath at least 1 minute. Then, add other components according to the table.

Table 2. Reaction Setup - First strand cDNA buffer (Mixture B) per reaction (This Master Mix can be prepared before or during the denaturing step)	
Components	Volume
5X RT Buffer (DTT/dNTPs)	4 μ l
DEPC-Treated H ₂ O	5 μ l
RTase/RI Enzyme Mix	1 μ l
Final volume	10 μl

3. Mix the following reaction mixture, then incubate at 25°C for 10 minutes and 37-50°C for 50 minutes.

Table 3. Reaction Setup	
Components	Volume
Mixture A	10 μ l
Mixture B	10 μ l
Final volume	20 μl

4. Incubate at 85°C for 5 minutes for Termination.
5. Store cDNA at -20°C or for immediate real-time PCR reaction

Recommended real-time PCR Condition

Table 4. Reaction Setup for real-time PCR	
Components	Volume
cDNA	2 μ l (100 fg -100 ng) *



Forward primer	50 – 400 nM**
Reverse primer	50 – 400 nM**
2X qPCR Master Mix	10 µl
H ₂ O	to a final volume of 20 µl
Final volume	20 µl

* The volume of cDNA should not be more than 10% of the total qPCR reaction volume.

** The PCR primer concentration for an optimal qPCR reaction may vary based on the properties and the template of the primers.

Recommended real-time PCR Program

Step	Temperature	Time	Cycles
Template denaturation and enzyme activation	95°C	10 mins [#]	1
Denaturation	95°C	15 sec	40
Annealing and Extension	60°C	60 sec	
Melting curve analysis	Refer to the instrument manual		

[#] 10 minutes is recommended for the first step to completely denature the DNA and activate the enzyme.

Recommended PCR Condition

Components	Volume
cDNA	2-10 µl
10X <i>Taq</i> Buffer	5 µl
Forward primer	0.1 – 0.5 µM
Reverse primer	0.1 – 0.5 µM
dNTPs	0.2 mM each
<i>Taq</i> DNA polymerase	0.25 µl (1.25 units)
H ₂ O	to a final volume of 50 µl
Final volume	50 µl



Recommended PCR Program

Step	Temperature	Time
Initial denaturation	94°C	2 mins
Denaturation	94°C	30 sec
Annealing	50-68°C ##	30 sec
Extension	72°C	30 sec/kb
Final Extension	72°C	1 min

} 25-40 cycles

The optimal PCR conditions differ based on the thermodynamic properties of the primers.

© **Revision History** ©

Description	Version	Date
Initial Release	S2GNM03j30002_Protocol_V1	Aug 2023