



10mM dNTPs Mixtures (2.5 mM each)

(Cat. No.: S2GNM03g30001)

Description:

The 10mM dNTPs Mixtures are supplied as ready-to-use solutions, including 2.5 mM each of dATP, dCTP, dGTP, and dTTP. This product saves researchers time in manual mixing and minimizes the possibility of contamination.

The deoxynucleotides are suitable for experimental applications in molecular biology, such as PCR, DNA sequencing, cDNA synthesis, etc.

Kit Contents:

Contents	S2GNM03g30001
10mM dNTPs Mixtures (2.5 mM each)	1mL

Storage:

Stable for 24 months at -20°C.

© Revision History ©

Description	Version	Date
Initial Release	S2GNM03g30001_Protocol_V1	June 2023