



40mM dNTPs Mixtures (10 mM each)

(Cat. No.: S2GNM03g30002)

Description:

The 40mM dNTPs Mixtures are supplied as ready-to-use solutions, including 10 mM each of dATP, dCTP, dGTP, and dTTP. This product saves researchers time in manual mixing and minimizes the possibility of contamination.

The deoxynucleotides are suitable for experimental applications in molecular biology, such as PCR, DNA sequencing, cDNA synthesis, etc.

Kit Contents:

Contents	S2GNM03g30002
40mM dNTPs Mixtures (10 mM each)	1mL

Storage:

Stable for 24 months at -20°C.

© Revision History ©

Description	Version	Date
Initial Release	S2GNM03g30002_Protocol_V1	June 2023