

**100mM dNTP Set (100mM each)****(Cat. No.: S2GNM03g30003)****Description:**

The 100 mM dNTP Set includes 100 mM each of dATP, dCTP, dGTP, and dTTP. Each vial contains 500 µl in purified water (pH 7.5).

The deoxynucleotides are suitable for experimental applications in molecular biology, such as PCR, DNA sequencing, cDNA synthesis, etc.

Kit Contents:

Contents	S2GNM03g30003
100mM dNTP Set	4X 0.5mL (100mM each)

Storage:

Stable for 24 months at -20°C.

◎ **Revision History** ◎

Description	Version	Date
Initial Release	S2GNM03g30003_Protocol_V1	June 2023